



Huachenbio -- Your Professional Choice

- ✦ Suitable for Primary Cell Isolation and Subculture
- ✦ Compatible with a variety of mesenchymal stem cells, such as umbilical cord-, adipose-, bone marrow-, amniotic membrane-, hair follicle-, and dental pulp-derived mesenchymal stem cells
- ✦ Serum-free, free of any animal-derived components, antibiotic-free, stable in performance, and minimal batch-to-batch variation
- ✦ High cell expansion rate, with a single passage expansion fold of over 20x
- ✦ Cell yield per T175 flask: $>2 \times 10^7$ cells;
Cell yield per 10-layer cell factory: $8-10 \times 10^8$ cells
- ✦ Cell diameter: 14-15 μm , smaller than that of similar products on the market
- ✦ GMP level, prepared with water for injection (WFI), endotoxin $< 0.1 \text{ EU/ml}$
- ✦ Independent R&D and production system, stable supply, and high cost-effectiveness



novastem-MSC[®]

间充质干细胞无血清培养基

Serum-Free Medium For Mesenchymal Stem Cell

苏州华辰生物科技有限公司

Suzhou Huachen Biotechnology Co., Ltd

Tel: 400-965-9800

Website: www.huachenbiotech.com

Address: Unit A3-504-1, Creative Industry Park, No. 328 Xinghu Street, Suzhou Industrial Park



WeChat
Official Account



WeChat Official
Subscription Account



WeChat ID
for Inquiries

COMPANY PROFILE

Suzhou Huachen Biotechnology Co., Ltd.

Founded by technical experts with years of in-depth experience in the Cell and Gene Therapy (CGT) field, Suzhou Huachen Biotechnology Co., Ltd. is a one-stop provider of full-process solutions for cell and gene therapy. It focuses on the R&D and sales of high-end reagents in the CGT sector, offering customers high-end reagents that can meet the needs of domestic substitution, as well as integrated process development solutions for cell therapy.

The company has obtained ISO 9001 Quality Management System certification, and some of its products have received Drug Master File (DMF) filing numbers from the U.S. Food and Drug Administration (FDA). It is recognized as a Leading Talent Enterprise in Suzhou Industrial Park and serves as a governing unit of the China Association for Medical Biotechnology.

It has now developed a series of product systems, including:

Hyper Series:

hyperClone NK Activation/Expansion Kit, hyperCryo Cell Cryopreservation Medium;

Nova Series:

novaStem-MSC Mesenchymal Stem Cell (MSC) Medium, novaNK-20 NK Cell Medium, novaT-15 T Cell Medium;

Star Series:

3D StarPore Microcarriers, StarPore Microcarrier-Specific Lysis Buffer, starMedium Chemically Defined Medium

PRODUCT INTRODUCTION

novastem-MSC Serum-Free Medium For Mesenchymal Stem Cell

It can be used for the primary cell isolation and subculture of mesenchymal stem cells (MSCs) from multiple sources, such as umbilical cord (UCMSCs), adipose tissue (ADSCs), bone marrow (BMSCs), amniotic membrane (AMSCs), hair follicles (HFSCs), dental pulp (DPSCs), etc., while maintaining the cells' pluripotent differentiation potential.

It is serum-free, free of any animal-derived components, and antibiotic-free, with stable performance and minimal batch-to-batch variation.

It features a high cell expansion rate, with the cell count at the same passage being significantly higher than that of similar products on the market.

It has an independent R&D and production system, ensuring stable supply and high cost-effectiveness.

It supports clinical-grade / pharmaceutical-grade cell culture.

USAGE INSTRUCTIONS >>>

01

- It is recommended that the supplement be thawed at 37°C.

Tips: Pay attention to the thawing time; it is optimal to leave a few small ice crystals when thawing is nearly complete. For small-volume use, aliquot and store the supplement in frozen state.

02

- Add it to the medium and mix thoroughly until uniform.


03

- Prepare a complete medium

It is recommended that the complete medium be prepared right before use and used up within one month.

If the culture system is small, it is advisable to aliquot the supplement into smaller portions for frozen storage based on actual usage. Prepare the complete medium in proportion when needed to avoid repeated freeze-thaw cycles.

Seeding Density



6000-7000 PCS/cm²



7000-8000 PCS/cm²



8000-9000 PCS/cm²



9000-10000 PCS/cm²

P1 P7

P8 P10

P11 P13

P14 and above

Cell Passage Time


It is generally around 3 days (72 hours). The proliferation rate of different human Mesenchymal Stem Cells (hMSCs) varies. It is recommended to determine the exact passage time based on cell confluency; passage is advised when cell confluency reaches approximately 80-90%. Excessive cell confluency (>95%) will affect subsequent cell growth.

Cell Morphology

The cells are spindle-shaped and exhibit a swirling (or vortex-like) growth pattern.


CULTURE OF MESENCHYMAL STEM CELLS

Primary Culture




- Isolate the Wharton's jelly tissue from the umbilical cord and perform primary culture using the tissue block adherence method. Inoculate 0.5g of Wharton's jelly into each T75 culture flask. After culturing in a 37°C cell incubator for 1-2 hours, add 6-8ml of cell culture medium. Supplement with 5ml of medium on days 2-3. Change the medium on days 7-10 by discarding the supernatant and adding 10-15ml of fresh medium. Observe the cell expansion on days 12-14 and proceed with digestion and passage.
- The earliest time for cells to migrate out of the tissue blocks is 5-7 days, and a small number of scattered cells can be observed under the microscope.
- The tissue block adherence efficiency is high, with a large number of tissue blocks adhering to the culture surface.
- A greater number of cells can be harvested.

Subculture



- Compared with similar products on the market, it has higher expansion efficiency. The cells at high passages age more slowly and can support culture for more than 20 passages.

Culture Performance



Passage	Seeding Density (cells/cm²)	Time	Culture Vessels	Harvested Cell Count (cells/flask)	Expansion Fold	Total Cell Number	Total Expansion Fold
P0	—	10–14天	T75	1.00E+06	—	2.00E+07	—
P1	7000	72h	T175	2.28E+07	18.6	3.72E+08	18.6
P2	7000	72h	T175	2.50E+07	20.4	7.59E+09	379.44
P3	7000	72h	T175	2.76E+07	22.5	1.71E+11	8537.4
P4	7000	72h	T175	2.79E+07	22.8	3.89E+12	194652.72
P5	7000	72h	T175	2.36E+07	19.3	7.51E+13	3756797.496
P6	7000	72h	T175	2.28E+07	18.1	1.36E+15	67998034.68
P7	7000	72h	T175	2.00E+07	16.3	2.22E+16	1108367965
P8	8000	72h	T175	2.10E+07	15.0	3.33E+17	16625519479
P9	8000	72h	T175	1.93E+07	13.8	4.59E+18	2.29432E+11
P10	8000	72h	T175	1.74E+07	12.4	5.69E+19	2.84496E+12
P11	9000	72h	T175	1.61E+07	10.2	5.80E+20	2.90186E+13
P12	9000	72h	T175	1.54E+07	9.8	5.69E+21	2.84382E+14
P13	9000	72h	T175	1.28E+07	8.1	4.61E+22	2.30349E+15
P14	10000	72h	T175	1.14E+07	6.5	2.99E+23	1.49727E+16
P15	10000	72h	T175	1.02E+07	5.8	1.74E+24	8.68418E+16

- The average expansion fold of the Novastem-MSC system from P1 to P5 is 20.7.
- Counted on average from 10 umbilical cords: for example, if approximately 10g of Wharton's jelly tissue is isolated, 2.0×10^7 cells can be obtained at Passage 0 (P0), and theoretically 1.71×10^{11} cells can be obtained at Passage 3 (P3).
- When prepared into products at the conventional dosage of 5×10^7 cells per unit, it can yield 3,420 units in total.

SUMMARY

Cell Passage Time



It is generally around 3 days (72 hours). Different hMSCs vary in their growth rates. It is recommended to determine the exact passage timing based on cell confluency; passage is advised when cell confluency reaches approximately 80%-90%. Excessively high cell confluency (>95%) will affect subsequent cell growth.

The medium volume is 0.2 mL/cm². No coating is required, and no medium change is needed for 3 consecutive days of culture.

When switching from other culture systems to novastem-MSD, the initial cell expansion fold may be relatively low. It is recommended to resuscitate and seed cells with a 1:1 mixture of the original medium and novastem-MSD. After 1 generation of culture, the system can be completely switched to novastem-MSD.

Cell Digestion



It is recommended to use a relatively gentle trypsin, such as Gibco TrypLE™ Express Enzyme (recombinant, 1X). Dilute it with PBS at a 1:1 ratio to prepare a working solution of trypsin before use. The digestion time is 3-5 minutes. After 1-2 minutes of digestion, gently tap the culture vessel; when most cells are observed to detach, terminate the digestion with medium.


Cell Cryopreservation



It is recommended that the final concentration of DMSO be 5%. Before cell digestion, you can first prepare a cryopreservation solution with a volume ratio of DMSO to complete medium of 1:9, and store it pre-cooled at 2-8°C. After the cells to be cryopreserved are mixed with complete medium, slowly add an equal volume of the pre-prepared cryopreservation solution, mix well, and then aliquot for cryopreservation.


QUALITY CONTROL ANALYSIS

Product Certificate of Analysis



After passing the release inspection, a corresponding Certificate of Analysis (COA) shall be issued. The main test items include characteristics, physical and chemical properties, safety indicators, and functional indicators.


Detection of Cell Surface Markers



Detection of Surface Markers in Cultured hMSCs

The positive rates of CD73, CD90, CD105, and HLA-ABC are >95.0%;
The positive rates of CD14, CD19, CD34, CD45, and HLA-DRDPDQ are ≤2.0%.

Cellular Multilineage Differentiation Potential



Three-Lineage Differentiation Assays:


A
Osteogenic differentiation
(Alizarin Red staining)

B
Adipogenic differentiation
(Oil Red O staining)


C
Chondrogenic differentiation
(Alcian Blue staining)

Detection of Cellular Immunomodulatory Capacity


CFSE Assay for Lymphocyte Proliferation




Inhibition of TNF-α Secretion




Upregulation Fold of Treg



Downregulation Ratio of Th1





Downregulation Ratio of Th17




The cultured human Mesenchymal Stem Cells (hMSCs) exhibit the following immuno-modulatory capabilities:

- Inhibiting lymphocyte proliferation;
- Promoting the upregulation of Regulatory T cells (Treg);
- Promoting the downregulation of T helper cell type 1 (Th1) and T helper cell type 17 (Th17);
- Inhibiting the expression of Tumor Necrosis Factor-α (TNF-α).

Detection of Cell Migration Ability
Detection of Angiogenic Ability




Karyotype Analysis of Cell Chromosomes



Karyotype Analysis Results of Continuously Cultured hMSCs:

No numerical or structural abnormalities of chromosomes were observed.

In Vitro Tumorigenicity Assay of Cells



A.P3 B.P8 C.P13

Results of In Vitro Tumorigenicity Assay for Continuously Cultured hMSCs:

Human Mesenchymal Stem Cells (hMSCs) failed to form colonies in the soft agar colony formation assay, indicating that they do not possess tumorigenic potential in vitro.

Ordering Information

Kit Name	Kit Cat#	Product Components	Product Cat#	Specs	Qty
Serum-Free Medium For Mesenchymal Stem Cell Kit	HC-UC02R5-kit	novaStem-MSCTM Serum-Free Medium For Mesenchymal Stem Cell(with Phenol Red)	HC-UC02R	500ml/bottle	1
		novaStem-MSCTM Basal Medium supplement	HC-UC02S	25ml/bottle	1
	HC-UC02F5-kit	novaStem-MSCTM Serum-Free Medium For Mesenchymal Stem Cell(without Phenol Red)	HC-UC02F	500ml/bottle	1
		novaStem-MSCTM Basal Medium supplement	HC-UC02S	25ml/bottle	1